Spontaneous electroretinogram signals to screen people with early risk factors for diabetic retinopathy

Ramsés Noguez Imm<sup>1\*</sup>, Julio Muñoz-Benitez<sup>2\*</sup>, Diego Medina<sup>2</sup>, Everardo Barcenas<sup>2</sup>, Guillermo Molero-Castillo<sup>2</sup>, Pamela Reyes-Ortega<sup>1</sup>, Jorge Armando Hughes-Cano<sup>1</sup>, Leticia Medrano-Gracia<sup>4</sup>, Manuel Miranda-Anaya<sup>6</sup>, Gerardo Rojas-Piloni<sup>1</sup>, Hugo Quiroz-Mercado<sup>7</sup>, Luis Fernando Hernández-Zimbrón<sup>7,8</sup>, Elisa Denisse Fajardo-Cruz<sup>7</sup>, Ezequiel Ferreyra-Severo<sup>7</sup>, Renata García-Franco<sup>9</sup>, Ellery López-Star<sup>10</sup>, Marlon García-Roa<sup>10</sup>, Van Charles Lansingh<sup>10</sup>, and Stéphanie C. Thébault<sup>1</sup>.

Corresponding author: Stéphanie C. Thébault, Ph.D. Email: <u>sthebault@comunidad.unam.mx</u>; \* Ramsés Noguez Imm and Julio Muñoz-Benitez contributed equally to this work.

<sup>1</sup> Instituto de Neurobiología y <sup>6</sup> Unidad Multidisciplinaria de Docencia e Investigación-Facultad de Ciencias; Universidad Nacional Autónoma de México (UNAM), Campus UNAM-Juriquilla, 76230 Querétaro, Mexico

<sup>2</sup> Facultad de Ingeniería; Universidad Nacional Autónoma de México (UNAM), Ciudad Universitaria, Ciudad de México, Mexico

<sup>4</sup> IIMAS; Nacional Autónoma de México (UNAM), CDMX, Mexico

<sup>7</sup> Research Department, Asociación Para Evitar la Ceguera, Mexico city, Mexico

<sup>8</sup> Clínica de Salud Visual, Escuela Nacional de Estudios Superiores, Unidad León, Universidad Nacional Autonóma de México (UNAM), León, Guanajuato, Mexico.

<sup>9</sup> Instituto de la Retina del Bajío (INDEREB), Prolongación Constituyentes 302 (Consultorios 410 y 411, torre 3, Hospital San José), El jacal, 76187, Santiago de Querétaro, Querétaro, Mexico

<sup>10</sup> Instituto Mexicano de Oftalmología (IMO), I.A.P., Circuito Exterior Estadio Corregidora Sn, Centro Sur, 76090 Santiago de Querétaro, Querétaro, Mexico

NOTE: This preprint reports new research that has not been certified by peer review and should not be used to guide clinical practice.

#### Abstract

Given the ever-increasing prevalence of type 2 diabetes and obesity, the pressure on global healthcare is expected to be colossal, especially in terms of blindness. Electroretinogram (ERG) has long been perceived as a first-use technique for diagnosing eye diseases, but existing methods are insufficient to screen early risk factors of diabetic retinopathy (DR). Here, we introduce non-evoked ERG as a simple, fast modality to record spontaneous activity, from which we developed a single random forest-based model that predicts disease cases in rodent models of obesity and in people with overweight, obesity, and metabolic syndrome. Classification performance was validated using a dataset from an independent eye center. Our algorithm can be coupled with different ERG sensors, including ones working with portative, non-mydriatic devices. Principal component and discriminant analysis suggest slow spontaneous ERG frequencies as main discriminators for our predictive model. Our study will facilitate the implementation of interventions for the prevention of overweight and obesity by providing a robust, quantitative, and non-invasive identification and follow-up approach, which should ultimately reduce DR incidence.

**Key words:** ERG, spontaneous activity, basal activity, resting-state activity, intrinsic activity, predictive biomarker, subclinical diabetic retinopathy, machine learning, random forest model

#### Introduction

From a public health standpoint, timely detection of diabetic eye disease is one of the most cost-effective health procedures available to avoid the burden of vision loss <sup>1</sup>. Huge efforts are being made to screen for sight-threatening diabetic retinopathy (DR) <sup>1,2</sup>, but the best intervention remains to prevent the onset of diabetes <sup>1</sup>.

DR affects nearly 100 million people worldwide <sup>3</sup>, corresponding to approximately one-third of all people with diabetes. Type 2 diabetes accounts for 85 %–95 % of diabetes cases and it can be avoided often by healthy life choices <sup>4</sup>. Strategies to promote adherence to such choices would benefit from reliable, large-scale screening methods that identify and allow follow-up of people at risk of type 2 diabetes, i.e., those with overweight, obesity or metabolic syndrome, but still without diabetes. Ultimately, predicting risk factors of developing type 2 diabetes should mitigate the risk of developing complications like DR. The screening method to be developed must overcome the limitations of the current gold standard tests, particularly the A<sub>1</sub>C <sup>5</sup> and the 2-hour post-challenge glucose tests <sup>6</sup> that both require fasting, and it should be based on a reliable biomarker of the silent installation of type 2 diabetes.

An appealing candidate is spontaneous retinal oscillations. Spontaneous neural oscillations have proven being biomarkers for neurodegenerative diseases <sup>7</sup> and diabetes is one of them <sup>8</sup>. Spontaneous brain activity changes with diabetes <sup>9</sup>, but also with obesity <sup>10–14</sup>. Notably, retinal neurons can produce spontaneous activities <sup>15–27</sup> and neurodegenerative retinas showed altered patterns of spontaneous activity <sup>28–30</sup>, including in diabetes conditions <sup>31</sup>.

Here we propose a diagnostic prediction method for early risk factors of DR based on the non-invasive recording of spontaneous retinal oscillations using a simple, yet meaningful non-evoked electroretinogram (ERG) protocol. ERG is the only clinical objective test

recommended by the International Society for Clinical Electrophysiology of Vision to stage early DR<sup>32</sup>, before vascular changes are apparent, and its clinical application has significantly improved with the commercial introduction of non-invasive, portative, and non-mydriatic ERG devices <sup>33</sup>. Nonetheless, ERG has always been based on the response of retinal cells to a flash of light and never under spontaneous conditions. We found that spontaneous ERG signals of 1-minute duration are differentially altered in rodent models of obesity and early diabetes, allowing their discrimination by a random forest-based prediction model. The model also predicts the evolution of the diabetes model and risk factors for DR in humans, including overweight, obesity, metabolic syndrome, and diabetes. Our algorithm can be coupled with spontaneous ERG signals from different sensors. Principal component and discriminant analysis revealed slow ERG frequencies as main discriminators for our predictive model. Together, our study shows that spontaneous ERG signals are intimately linked to systemic metabolic status and demonstrates their use to screen people for preventable stages of DR.

#### Results

# Predictive model for obesity, type 1 diabetes, and type 1 diabetes evolution based on spontaneous ERG oscillations in rodents

Spontaneous ERGs of high-fat diet-fed mice that are obese and insulin-resistant after 12 weeks (Supplementary Fig. 1A) showed three consistent peaks in the 0.1-10 Hz range (Fig. 1A) of similar power compared to lean mouse ERGs (Fig. 1B). The mid-low (1-1.8 Hz) peak frequency was reduced in obese mice, while the low (0.1-0.8 Hz) and mid (2-4 Hz) peak frequencies were not affected (Fig. 1B). Spontaneously obese insulin-resistant (Supplementary Fig. 1B) mice and streptozotocin-treated hyperglycemic (Supplementary Fig. 1C) rats exhibited low to mid-low frequency oscillatory activities (Fig. 1C, 1E, respectively) of similar power compared to their respective controls (Fig. 1D, 1F). The low (0.6-1 Hz) peak frequency (Fig. 1D) and the very low (0.2-0.6 Hz) and low-to-mid (0.6-2.5 Hz) peak frequencies (Fig. 1F) were reduced in ERGs of spontaneously obese mice and streptozotocin-treated rat, respectively. In contrast, oscillatory potentials (OPs), considered as the most precociously altered ERG parameter in diabetes <sup>34</sup>, remain unchanged in high-fat diet-fed mice (Supplementary Fig. 1D, 1E). We next tested the potential predictive content of these signals measured under the spontaneous modality. A single predictive model based on random forest algorithm was developed using the power spectra of spontaneous ERG oscillations in the activity range between 0.1-10 Hz (see Methods for more details). The area under the receiver operating characteristic (ROC) curve (AUC-ROC) values are 0.804, 0.875, and 0.906 for the high-fat diet-induced obesity mouse model (Fig. 1G), the spontaneously obese mouse model (Fig. 1H), and streptozotocin-induced type 1 diabetes rat model (Fig. 11), respectively. Additional metrics including accuracy and precision indicate that our model predicts correctly in 80 to 87.5 % of cases and is sure of its prediction in 80 to 84.6 % of

cases (Fig. 1G-1I). The specificity of this model is good (>0.75), meaning that it misses few non-diseased cases, and it is extremely sensitive in predicting both obesity and type 1 diabetes cases (Fig. 1H, 1J).

Predicting a disease condition is useful for diagnosis at a given disease stage, but diabetes and its complications evolve <sup>31,35</sup>. Progressive diabetic retinal neurodegeneration has been reported between 2 and 32 weeks of age in rodent models of diabetes <sup>36</sup>. We therefore asked if our model could discriminate the disease cases from control ones over time by using the power spectra of rat ERGs after 4, 6, 8 or 12 weeks of streptozotocin treatment. The performance of this multiclass classification is visualized in a confusion matrix (**Fig. 1J**). Our model correctly identified the week classes in 75 % of cases (**Fig. 1J**). These data show that the major earliest risks of DR, i.e., obesity and type 1 diabetes, can be predicted by spontaneous ERG oscillations in rodents.

#### Exceeding expert performance to predict preventable risks of diabetic retinopathy

To assess the clinical relevance of our approach, we created a human database of spontaneous ERGs (**Fig. 2A**) from 80 metabolically healthy adult subjects (54 eyes), 40 patients with overweight (24 eyes), 14 patients with obesity (8 eyes), 66 patients with metabolic syndrome but no diabetes (MetS, 56 eyes), and 68 patients with diabetes but no DR (44 eyes) (**Table 1**), who were referred to both IMO and INDEREB ophthalmology clinics for general eye check-up. Participants with MetS were older than controls (median 60.45 years vs 41.49 years, p<0.001) and patients with diabetes were older than all other groups (median 65.86 years vs 46.05 years, p<0.001), and more likely to be women (79.41 % vs 20.58 %, P < 0.0001) (**Table 1**). A ground-truth label for the presence of metabolic health, overweight, obesity, MetS, diabetes mellitus, and no DR was established by experts using internationally

accepted criteria, as referenced in the Methods. Importantly, OP parameters of the overweight, obese, and MetS groups remain unchanged compared with the control group and overall, they also showed no difference in the diabetes without DR group (Supplementary Fig. 2A-C). The DR score <sup>37</sup> is similar between the control, overweight, obese, and MetS group, and is reduced patients with diabetes but no DR (Supplementary Fig. 2D). We further showed that the spontaneous ERG oscillation component was neither correlated with OP parameters (Supplementary Fig. 2E-F) nor DR score (Supplementary Fig. 2G). Spontaneous activity is seen between 0.3 and 40 Hz (Fig. 2B), with main peaks in the 0.3-2, 10-20, and 20-40 Hz bands (Fig. 2B, insets). The large variation in our dataset that we assumed was confirmed by the absence of statistical differences in AUC and peak frequency between all groups in the previously mentioned bands (Fig. 2B). We therefore pooled spontaneous ERG-power spectra from non-healthy participants as one disease group (Fig. **2**C) and found an increased peak frequency in the 20-40 Hz band in the disease group (Fig. 2D). No significant difference was observed in the AUC or peak frequency in the 0.3-2 and 10-20 Hz bands (Fig. 2D). Next, the power spectra from 986 aleatory fragments of oneminute duration extracted from the 186 ERGs registered from a total of 268 patients (186 eyes), were randomized across training (80 %) and test (20 %) sets; diseases cases were defined as all cases except metabolically healthy and no DR, and the same random forestbased model used in animal models was applied in our database of human ERG power spectra.

The predictive model performed good (AUC-ROC = 0.761, **Fig. 2E**) and is accurate in 67.3 % of the cases (**Fig. 2F**). The prediction precision is poor and the model shows 85.4 % sensitivity and a specificity of 61.5 % (**Fig. 2F**). We then aimed to assess the accuracy of our model using an external cohort. To this end, we created a validation database based on data

from patients enrolled in an important eye care setting in Mexico City (**Table 2**). Using this validation dataset of 0.3-40 Hz power spectra (126 aleatory fragments of one-minute duration extracted from the 38 ERGs registered from a total of 39 patients (11 eyes from controls, 5 eyes from overweight, 3 eyes from obese, 3 eyes from MetS, and 3 eyes from diabetes but no DR), the predictive model performance was similar to the one of the original model (**Fig. 2G, H**). These data provide an external validation of our predictive model for early DR risks. Interestingly, when we fed our model with the power spectra from ERGs registered with one and same device (762 aleatory fragments of one-minute duration extracted from 173 ERGs registered from a total of 237 patients (**Table 3**, 173 eyes)), the predictive model performed greatly (AUC-ROC = 0.926, **Fig. 2J**) and is accurate in 86.6 % of the cases (**Fig. 2K**). The prediction precision is 89.8 %; the model shows 82.8 % sensitivity and a specificity of 90.4 % (**Fig. 2K**).

During the model implementation, we found that the random forest-based model outperformed other models based on other algorithms, including linear and radial support vector machines (**Supplementary Fig. 3A, B**) and deep learning (**not shown**). We also sought the best performance for the shortest ERG window to shorten test duration. The greatest prediction performance was established for spontaneous photopic ERGs of 60 s (**Supplementary Fig. 3C-E**).

Case labeling relies on expert's knowledge. When we compared our model predictions with the predictions made by medical experts, our model underperformed in terms of precision (0.361) but was consistent in recall (0.867) (**Fig. 2I**). This, when experts only access full eye examinations (**Fig. 2I**). Our model performance slightly improved (precision 0.444) in case experts had access to reference tests to diagnose potential diabetes (**Fig. 2I**). Our model was close to the judgements of experts in ophthalmology if, in addition to eye examination, they

had access to laboratory (glycemia, HbA1c levels, HOMA-I, and lipid profile) and nonlaboratory (blood pressure, weight and waist circumference, and body mass index (BMI)) tests (**Fig. 2I**). In the three contexts, our model showed a good recall, referring to few false negatives –cases our model would classify as healthy that are not– (**Fig. 2I**). The fact that our model recognized disease cases with less information than experts prompted us to look closely at cases of errors in our model in all three contexts (only complete eye exam, the latter plus reference blood test; the latter two plus specific tests for metabolic profile). It appeared that error cases are false positives (low precision), indicating that our model "incorrectly" labeled cases as positive (or diseased) that were labeled as negative by experts. But are these negative labels from experts truly negative? Experts classified them as "with no retinopathy" based on eye examination, while deeper blood tests revealed that these "negative cases" present risk factors for DR. These data demonstrate that our model can predict early risk factors for DR with less information than experts, outperforming them.

# Statistical explanation for the spontaneous ERG-based model of prediction for preventable risks of DR

When we implemented a multiclass model encompassing all risk stages for diabetes and DR, classification performance reached 66 % (**Fig. 3A**). The largest errors of our model consisted in classifying obesity into more advanced states and some overweight cases into controls (**Fig. 3A**).

To gain insights into the discrimination between groups, we performed explanatory statistical analyses. We considered the variables with higher variability (detailed in Methods), including the normalized power of the slowest oscillations (0.3-0.8 Hz), the AUC of the 0.3-40 Hz band, and the peak frequency of the 20-40 Hz band. Measurements of the considered

plasma and body variables significantly differed between control and disease groups (Table 1 and 2), no significant change in the power of slow frequency spontaneous oscillations of the ERG was detected, the peak frequency in the 20-40 Hz band increased in the disease group (Fig. 2B, D). The initial evaluation of the utility of these variables in separating control and disease groups was accomplished by examining their natural partitioning using Principal Component Analysis (PCA). When only using variables related to spontaneous ERG oscillations, no separation of all groups was observed, but control cases tend to cluster (Fig. **3B**). PC1 and PC2 explained 69 % of the data set variance, powers of 0.6-0.7 Hz correlate with the AUC of the 0.3-2 and 0.3-40 Hz bands and had the largest contributions to PC1 (Fig. 3C). The separation of the control group could be better appreciated when spontaneous ERG oscillation-related variables were combined with body variables like triglyceridemia, total cholesterol levels, fasting blood glucose, very low-, low-, and high-density lipoprotein cholesterol, systolic and diastolic blood pressure, age, body weight, hip and waist circumferences, and insulinemia. PC1 and PC2 explained 47 % of the data set variance (Fig. **3D**). According to our PCA results, the normalized power of the slowest oscillations (from 0.4 to 0.7 Hz) was the most important in separating disease cases from controls as they had the largest contributions to PC1 (Fig. 3D). Plasma triglycerides, total cholesterol, low- and high-density lipoprotein cholesterol, fasting blood glucose, insulin, systolic blood pressure, age, and hip circumference, were the least important in this regard (Fig. 3D). Based on visual inspection of the PC1/PC2 biplot (Fig. 3E), we refined our PCA analysis to variables related to spontaneous ERG oscillations (from 0.7 to 0.8 Hz, the peak frequency in the 20-40 Hz band, and the AUC of the 0.3-40 Hz band) that were non-orthogonal with body-related variables (Age, fasting blood glucose levels, and total cholesterol), but PCA did not show a better separation of disease and control cases (Fig. 3F). We finally used linear discriminant

analysis since it is geared towards the discrimination between user-defined groups <sup>40</sup>. As shown in the confusion matrix of the best performances we found for discriminating control (76 %), risk factors for diabetes (78 %), and diabetes (60 %) cases (**Fig. 3G**), this relative success required the combination of body metrics in addition to spontaneous ERG oscillation-related variables.

#### Discussion

The clinical community is experiencing an explosion of machine learning-guided (tele)diagnostics, particularly in the context of detecting treatment-requiring DR, and the positive impact on care delivery is becoming evident <sup>41,42</sup>. Risk-scoring algorithms for undiagnosed diabetes are also available <sup>43</sup>, but there is no such approach for the main risk factors of type 2 diabetes, i.e., overweight, obesity, and MetS. The purpose of screening such conditions resides in well-established data about their preventability with adopting healthy lifestyle behaviours <sup>44</sup> and because people of normal weight are not always metabolically healthy <sup>45</sup> and at the opposite, the obesity phenotype can associate with no or little evidence of metabolic dysfunction <sup>46</sup>. We have exploited the potential of non-invasive, quantitative, and objective ERG by introducing a new and simple modality to record spontaneous activity of the retina and combined it with a supervised machine learning-based model to predict early risk factors of DR. Our data support the robustness of the random forest system in the screening of early risk factors for DR and imply its translation into clinical use thanks to automated platforms like the one we created (http://deepretinopathydx.inb.unam.mx/). Our findings also highlight a very early impact of systemic metabolic changes in spontaneous signals from the central nervous system and add to the growing list of evidence that show retinal neurodegeneration as an early event in DR pathogenesis <sup>47</sup>.

#### Strengths and interpretation of spontaneous ERG signals' predictive content

We found that the spontaneous oscillations detected by non-evoked ERG are a quantifiable biological parameter that is modified under conditions related to an excess of body weight. These alterations were not accompanied by changes in OPs in high-fat diet-fed mice and in people with overweight, obesity, and MetS. These observations are consistent with previous

studies in animal models <sup>48,49</sup>, but contrast with data showing reduced OP amplitude in highfat diet-fed mice <sup>50</sup> or in ob/ob mice <sup>51</sup>. These discrepancies may be attributed to the different induction mechanisms of obesity, known to result in models with their own characteristics <sup>52,53</sup>. Particularly, our high-fat diet is more enriched in lipids than carbohydrates <sup>50</sup>. While the ERG B-wave is affected in patients with obesity <sup>54,55</sup>, there is no previous report of the effect of obesity on OPs in humans. In our group of patients with diabetes but no DR, the amplitude and implicit time of the OPs tend to decrease and increase, respectively, as previously reported <sup>56</sup>. In general, these data agree with a view that the spontaneous activity of the retina may replace OPs by becoming <sup>34</sup>, the most precociously altered functional parameter in prediabetic stages. In further agreement with this is our finding that the spontaneous ERG oscillation-related component score does not correlate with either OP parameters nor DR score that reflects the risk of requiring intervention within 3 years <sup>37</sup>.

The OPs are thought to reflect the function of the inner retina and are sensitive to changes in retinal circulation <sup>56</sup>. In this sense, one may wonder about the nature of spontaneous oscillations measured by non-evoked ERG. If ERG is classically conceived as the summation of local synaptic and intrinsic activities of retinal cells, the low amplitude and slow frequency oscillations recorded in the non-evoked mode may be contaminated by sources other than the retina, e.g., cardiorespiratory system and brain. In humans, with our setup setting, breathing is too slow (~0.2-0.3 Hz) <sup>57</sup> to be detected by ERG, as is eye movement (0.26 Hz) <sup>58</sup>, but cardiac (0.9-1.4 Hz) <sup>59</sup> and brain activities may be. Also, in stably anesthetized mice, respiration is in the 0.9-1.08 Hz range, while cardiac rhythm lies between 5 and 7.5 Hz. In ketamine-xylasine-anesthetized rats, eye movement shows two frequency components, 1 and 12.2 Hz, the former correlating with respiration frequency <sup>60</sup>. We cannot ascertain that the full range of spontaneous oscillations measured by non-evoked ERG originates in the retina,

but it is plausible that at least part of these signals come from the retina because in healthy adult retinas, several types of neurons have been reported to spontaneously oscillate, and this, within a range of frequencies from 0.7 to >10 Hz. In particular, spontaneous  $Ca^{2+}$ -dependent membrane oscillations have been recorded in bipolar cell axon terminals <sup>22,23,61</sup> and if the consequent pulsatile release of neurotransmitter drives rhythmic activity in post-synaptic neurons, including amacrine and ganglion cells<sup>21,62</sup>, different types of amacrine cells are also able to produce intrinsic oscillatory activity <sup>24,25</sup>. Low-amplitude oscillations have even been recorded in starbust amacrine cells<sup>21</sup>. The presence of these intrinsic oscillators in the inner retina <sup>21–25,61,62</sup>, the fact that retinal neurons can be electrically coupled <sup>63</sup>, and the large-scale network interactions happening in the retina<sup>20</sup> are likely to result in spontaneous fluctuations of the field potential <sup>64</sup> in this tissue. Consistent with the existence of bipolar/amacrine cell oscillators and our findings showing that spontaneous activity is altered in the high-fat dietinduced obesity, spontaneous obesity, and streptozotocin-induced type 1 diabetes models, inner retinal deficits have been detected at the onset of diabetes <sup>31,65,66</sup>. Overall, we found reduced peak frequencies of the low range activities (0.2-2.5 Hz), which is agreement with decreased inhibition in the early diabetic retina <sup>31,67-71</sup>. Also in favor of part of the spontaneous oscillations detected by the non-evoked ERG being produced in the retina is the recent finding that DR can be detected by machine learning processing of electrooculogram (EOG) signal —corresponding to the potential difference between cornea and retina—<sup>72</sup>.

Though studies of spontaneous retinal activity are sparse in the adult and its functional role is yet to be fully understood <sup>73</sup>, our results support that the 0.3-40 Hz activity is relevant for the disease process. Our data show that a single model is able to distinguish disease cases from control cases under disease-relevant conditions that range from rodent models to patients. This suggests that the predictive content of spontaneous ERG signals is conserved

in mammals and robust. Our PCA analysis adds to this view, because it showed that the normalized powers of the slowest oscillations (0.4-0.7 Hz) were the most important variables in separating disease cases from controls. Nevertheless, changes in spontaneous ERG oscillations of higher variability were insufficient to separate the groups of interest using PCA. Compared to the PCA that does the separation in a single step, the Random Forest model does it in several stages, that is, it first classifies the disease and control cases and then, in the disease category, classifies the four groups of interest, which surely contributes to its performance compared to that of PCA. There is also the algebraic limitation of PCA which cannot process more variables (823 in total) than the number of individuals (375, Table 1 and Table 2). We do not exclude the possibility that the frequency variables analyzed by PCA may explain the multigroup separation once the number of participants exceeds the number of variables. The statistical analyses we undertook for an explanation gave us some hints: the normalized power of the slowest oscillations (from 0.3 to 0.5 Hz) and the peak frequency of the 20-40 Hz band did not correlate with the power of the 0.6-0.8 Hz components or with the AUC of the 0.3-2 and 0.3-40 Hz bands. However, we do not have yet enough data to understand in details how our predictive model uses spontaneous ERG data and the fact that it surpassed so-called "black box" algorithms such as deep learning, does not yet represent an advantage in terms of intelligibility of our model <sup>74</sup>. If we cannot explain the specific loadings of each variable and we do not exclude the contribution of other variables, our findings support a very tight relationship between systemic metabolic changes and retinal function 75.

Our predictive diagnostic system meets the Food and Drug Administration criteria in terms of precision ( $\geq$  85 %), sensitivity ( $\geq$  85 %), and specificity ( $\geq$  82.5 %) for clinical validation of diagnostic tools in retina <sup>76,77</sup> and the validated ROC of 0.719 alsmot indicates excellent

clinical accuracy <sup>78</sup>. It performed as good as experts in identifying disease cases, but without the need for invasive study data. Our PCA data abound in this direction since additional variables to the spontaneous oscillations of the ERG such as VLDL, body weight, diastolic blood pressure and waist circumference, are necessary to better separate the control from disease groups. Importantly, accuracy is maintained when using signals from different sensors, which also accounts for the robustness of our system. It is well known that data from different sources penalize the performance of prediction models <sup>79</sup>. Moreover, the fact that spontaneous ERG protocols could be developed on three commercial sensors and that the predictive potential of these signals is shown here illustrates that our system is flexible and can be adapted to any ERG apparatus. Also, we found that predictions are better with photopic recordings, which is advantageous in clinical practice because they take less time than the dark-adapted ones. Based on our practical experience, we specify that 2 to 5 minutes of high signal/noise ratio recordings are enough to obtain the spontaneous 1-minute ERG sequences that, according to our data, are informative. Furthermore, the notion that the spontaneous oscillations detected by non-evoked ERG are the most precociously altered functional parameter in prediabetic stages is highly relevant in clinical practice because OPs can only be extracted from scotopic ERG <sup>80</sup>. Therefore, having access to such a parameter without the need for adaptation to darkness and light flash represents a definite advantage for the patient (no pupillary dilalation, faster examination <sup>80</sup>). For the above and because it considers existing resources in terms of ERG device, works with completely non-invasive, portatile and cost-effective ERG devices, and does not require ophthalmological experts, our system has easy applicability in clinical settings.

The clinically relevant information our system contributes to is to detect people with high likelihood of preventable risk factors for DR. Even though the Diabetes Prevention Program

recently reported that interventions that delay the development of type 2 diabetes in those at risk (overweight/obese with dysglycemia) do not reduce the subsequent prevalence of DR<sup>81</sup>, this concerns people who have progressed to diabetes. In cases where intensive lifestyle intervention and metformin managed to prevent progression to type 2 diabetes, DR does not happen <sup>81</sup>. In this regard, it is important to recognize that lifestyle interventions must be carried out with considerable involvement of clinicians and that acheiving awareness among patients and family relatives is challenging.

#### Limitations and future work

The advantages of our predictive diagnostic system —such as increased objectivity and efficiency in determining early risk factors for DR by the Random Forest system compared with health-care professionals, higher referral adherence from real-time point-of-care screening recommendations <sup>42</sup>, more efficient resource allocation towards prevention and treatments due to the Random Forest system offloading tasks from human graders, and reduction in the prevalence of RD in the mid to long term— are implied but unproven. Also, although translation of our system into clinics seems feasible, real-world testing and acceptation is in its infancy. We have provided external validation of our system, but future data from multiethnic populations should consolidate the validation of our model. In the same line, machine learning models have the disadvantage of dealing with the problem on which they were trained. If our findings in the streptozotocin-induced diabetes model suggest that our predictive model may be useful for disease follow-up, we plan to train our model with more severe grades of DR and other retinal disorders, step by step, possibly finishing with a multiclass prediction with acceptable performance in an unselected population. This will also

help clarify if altered spontaneous ERG oscillations represent neural signatures specific for each neurodegenerative retinal disease.

An additional crucial step is to compare our system with more gold standards of early DR, like the multifocal ERG implicit time <sup>82</sup>, and to determine whether the spontaneous 0.1-10 Hz oscillations in rodents and between 0.3-40 Hz in humans respond to a therapeutic intervention <sup>83</sup>. Showing altered peak frequencies and high degrees of prediction of these features only suggests that they are necessary for the early process of the diabetic eye desease <sup>84,85</sup>. Deciphering pathological mechanisms responsible for these alterations will be highly informative about DR etiology and may benefit the development of therapeutic options. Additional work is also needed to address the mechanisms that govern the slow spontaneous ERG signal.

#### **Concluding remarks**

Our ultimate goal is to provide a large-scale, easy, and affordable screening method that identifies asymptomatic patients in avoidable stages to enable personalized preventive action against the main early causes of DR. The state-of-the-art performance of our unique approach will likely contribute to improving the reputation of ERG <sup>86</sup>, putting it at the right place in the clinical scene, as a quick, easy to administer and interpret, and relevant tool for screening early DR risk factors.

#### Methods

#### **Ethics statement**

All animal experiments were approved by the Bioethics Committee of the Institute of Neurobiology (protocol #74) at UNAM (clave NOM-062-ZOO-1999) in accordance with the rules and regulations of the Society for Neuroscience: Policies on the Use of Animals and Humans in Neuroscience Research. Approval was obtained from the IMO and INDEREB Human Participants Ethics committee (reference: CEI/029-1/2015), the National Ethics Committee (reference: CONBIOÉTICA-09-CEI-006-20170306), the Research Committee at APEC (17 CI 09 003 142), and the Research Ethics Committee at ENES León (reference: CEI\_22\_06\_S21). Written informed consent was provided by all subjects. All procedures were conducted in accordance with the tenets of the Declaration of Helsinki.

#### Animals and models

C57BL/6 mice (male:female ratio = 1) between 6 and 8 weeks of age were obtained from the Institute of Neurobiology's animal house. Male Wistar (n = 20-22) adult (250-300 g) rats (n = 6-8) were used, as well as lean (n = 15) and spontaneous obese (n = 12) *Neotomodon alstoni* mice <sup>88</sup>. Animals were fed *ad libitum* and reared in normal cyclic light conditions (12 h light/dark cycle) with an ambient light level of ~400 lux. Plasma glucose concentrations were measured from a tail blood sample using a blood glucose analyzer (Accu-check active, GC model).

Six- to eight-week-old C57BL/6 mice were divided into two groups of 16 and fed a chow diet (5020, Lab Diets) containing 21 % of calories from fat or a high-fat diet containing 60 % calories from fat (D12492 Research Diets).

Diabetes was induced in Wistar rats by i.p. injection of streptozotocin (55 mg per kg body weight) in citrate buffer <sup>89</sup>. Control rat group received only citrate buffer injections. We confirmed diabetes by measuring blood glucose (>250 mg/dl in animals <sup>90</sup>) 24 h after streptozotocin injection. Bodyweight was measured as indicated and glycemic controls were always performed after a 6 h fasting <sup>91</sup>.

For in vivo ERG, animals were evaluated after a 12-h dark adaptation period.

#### Insulin and glucose tolerance tests

Glucose tolerance tests and insulin tolerance tests were performed on C57BL/6J after 12 weeks of control or high-fat diet and in 1-year-old lean and spontaneously obese *Neotomodon alstoni* mice, after 6-h fasting <sup>91</sup>.

The insulin tolerance test consisted in measuring glucose levels with a glucometer in tail vein blood samples obtained with a lancet needle before or 15, 30, 45, 60, and 90 minutes after an ip injection of 1 U/kg insulin (Humulin R; Eli Lilly). For the glucose tolerance test, mice were given glucose at a dose of 2 g/kg through a jugular vein catheter. Blood samples were then collected at 0, 15, 30, 45, 60, and 90 min after glucose administration to measure glycemia. Glucose profiles normalized to the initial glucose reading of each mouse were plotted for each group versus time of subsequent glucose determinations.

#### **Electroretinograms in animals**

Animals were anesthetized with 70 % ketamine and 30 % xylazine (1  $\mu$ l/g body weight, ip). Corneas were hydrated with hypromellose (5 mg/mL), and pupils were dilated with tropicamide-phenylephrine (50 mg/8 mg/mL). ERG responses were recorded with contact

lens silver electrodes <sup>92</sup> (3.0 mm diameter, Ocuscience) placed on each cornea. Reference electrodes were positioned subcutaneously exactly between the eyes.

The signal was amplified x100; the bandpass was set at 0.1 Hz to 1 kHz (AC-DC Differential Amplifier, Model 3000, A-M Systems) and acquisition frequency to 1 kHz (USB-6009, National Instruments). Spontaneous mesopic activity was measured for 10 minutes after a 20-minute dark adaptation period (0.1 lux) and then spontaneous photopic activity for 10 minutes after adaption to normal light (400 lux) for 10 min. At the end, light stimulation (0.7 ms flashes of 0.38 log cd.s/m<sup>2</sup>; MGS-2 white Mini-Ganzfeld Stimulator, LKC Technologies) was given to confirm retina function. If no classical evoked response was seen <sup>93</sup>, data were discarded.

#### Human data description

A total of 520 adult subjects aged between 30 and 80 years (mean:  $45.27 \pm 0.82$  years, 265 females), metabolically healthy or with overweight, obesity, MetS, or diabetes but no DR, were enrolled between February 26, 2015 and December 2019 and from September 2021 and June 17, 2022, in the Mexican Institute of Ophthalmology (IMO) of Querétaro (mean age:  $51.39 \pm 1.49$  years, 27 females), between November 11 and December 20, 2019 and from January 6 and May 26, 2022, in the Instituto de la Retina del Bajío (INDEREB) in Querétaro (mean age:  $32.98 \pm 2.07$  years, 24 females), between August 10, 2021 and March 20, 2022 in the Asociación Para Evitar la Ceguera (APEC) in Mexico city (mean age:  $45.77 \pm 1.20$  years, 119 females), and between August 10, 2021 and May 31, 2022 in the Clínica de Salud Visual (CSV) at ENES León UNAM in León (mean age:  $44.96 \pm 0.84$  years, 96 females). 375 (age mean:  $46.01 \pm 0.98$  years, 180 females) completed all tests required for the current study.

Subjects underwent an anamnesis and an initial optometric examination to ensure that they were eligible to participate. The exclusion criteria were ages outside 30 to 80 range, lens opacity, myopia greater than 6 diopters, glaucoma or other concomitant ophthalmologic disorders, ocular anomalies (e.g., surgery, trauma), recent use of laser or anti-angiogenic intravitreal administration, and cornea problems that disable ERG recordings.

Qualified medical personnel collected the anthropometric data in the morning (8 am to 9 am) after an overnight fast. Height was measured to the nearest 0.5 cm with a stadiometer (Seca 213; Seca). Bodyweight was measured with subjects wearing light clothing and without shoes to the nearest 0.1 kg on a mechanical column scale (Seca 700; Seca). Waist circumference was measured on undressed subjects at the midpoint between the lower margin of the last palpable rib and the top of the iliac crest while the subject was standing, after a moderate expiration, with a non-stretchable tape. Blood pressure was measured by using a mercury cuff sphygmomanometer after the study participant had been quietly seated for  $\geq 10$ min. Blood samples were taken from an intravenous catheter without stasis after an overnight fast of at least 8 h. Laboratory measurements that include fasting blood glucose, glycated hemoglobin (HbA1c), insulin, triglycerides (TG), low-density lipoprotein (LDL) cholesterol, very low-density lipoprotein (VLDL) cholesterol, high-density lipoprotein (HDL) cholesterol, total cholesterol (CT), and creatinine were performed at INTERMEDIC (Querétaro, Mexico) for IMO data, iml Laboratorio Médico (Querétaro, Mexico) for INDEREB data, Laboratorio clínico Jenner (Mexico City, Mexico) for APEC data, and Laboratorios Salud Digna (León, Mexico) for CSV data. The homeostasis model of assessment index (HOMA-I) was calculated using fasting insulin and glucose values <sup>94</sup>. All patients were classified according to the following criteria: normoglycemia (fasting glucose < 6.1 mmol/l) and diabetes (fasting glucose  $\ge 7.0 \text{ mmol/l}$ ), according to the

1999/2006 WHO criteria <sup>95,96</sup>. Normal weight was defined body mass index (BMI) between 18.5 and 24.99 kg/m<sup>2</sup>, overweight with a BMI between 25 and 29.99 kg/m<sup>2</sup>, and obesity with a BMI over 30 kg/m<sup>2</sup>. MetS was defined according to the International Diabetes Foundation criteria for MetS <sup>97</sup>, i.e. central obesity —BMI>30 kg/m<sup>2</sup> or waist circumference >= 94 cm in male and >= 80 cm in female— plus any two of the following four factors: raised TG, reduced HDL, raised blood pressure, raised fasting plasma glucose, raised HbA1c, raised plasma insulin, TG, LDL, VLDL, total cholesterol, creatinine, HOMA-I, or atherogenic index (calculated as log<sub>10</sub> (TG/HDL) <sup>98</sup>). Data from IMO, INDEREB, and APEC (*n* = 307) were used for model training and test, and data from CSV (*n* = 39) were used for external validation of the predictive diagnosis model. **Table 1** and **Table 2** contain an overview of patient demographics and biometrics for training/test and external validation of the predictive diagnosis model, respectively.

All subjects underwent a complete ophthalmologic examination including visual acuity testing using Snellen primer; anterior segment and crystalline status under microscopy and indirect ophthalmoscopy with a magnifying glass of 20 diopters; intraocular pressure by flattening tonometry (iCare TA01i); photographic study (7-field color photographs under pupil dilation, ZEISS camera, FM/NA, 60-degree images at IMO; ZEISS clarus® 500 Fundus Camera at INDEREB and CSV; and Visucam® 500 at APEC; and macular patterns, raster and macular thickness map by optical coherence tomography (OPTOVUE RTV-1000 equipment at IMO; Spectralis® Heidelberg Engineering at INDEREB and APEC; and CIRRUS HD-OCT 5000, Zeiss at CSV). Ophthalmic ERG tests were also performed, following ISCEV guidelines <sup>80</sup>. At IMO, the dark-adapted 0.01, 3.0, and 10, dark-adapted 3.0 oscillatory potentials, light-adapted 3.0, light-adapted 3.0 flicker, and multifocal ERG responses were measured in the order indicated using the MonElec2 (Metrovision, France;

from February 26, 2015 to December 20, 2017) or Retimax Advanced (CSO, Italy; from January 10, 2018 to June 17, 2022). At INDEREB, APEC, and CSV, the dark-adapted 0.01, 3.0, and 10, light-adapted 3.0, light-adapted 3.0 flicker ERG responses and the DR score <sup>37</sup> were recorded using a mydriasis-free ERG device (RETeval complete, LKC Technologies, USA). OP features (N1, P1, and N2 amplitudes and peak times, N1-P1 and P1-N2 ratios) were extracted from ISCEV dark-adapted 3.0 ERG protocol <sup>80</sup>. The DR score derived from flicker ERG and pupillography data correlate with ocular intervention for DR <sup>37</sup>. In addition, spontaneous ERG responses were measured in all patients using a custom protocol developed specifically for each ERG device (no flashlight, 0.3-1000 Hz band-pass filter with a 50 Hz notch, 1 kHz acquisition, and ×100,000 gain). At IMO, 30 min prior to ERG recording, one drop of tropicamide 1 % was instilled into each eye as a cycloplegic and contact lens electrodes were used posterior to corneal anesthesia with proparacaine hydrochloride eye drops. At INDEREB, APEC and CSV, skin electrodes on the lower eyelid were used <sup>99</sup>. All patients were adapted to mesopic conditions for 20 minutes prior dark-adapted ERG assessment and before ERG assessment under photopic conditions, patients were adapted to normal light (400 lux) for 10 min. The acquisition sequence was as follows: spontaneous dark-adapted ERG, evoked dark-adapted ERGs, spontaneous light-adapted ERG, and evoked light-adapted ERGs. Electrode impedance was monitored throughout the test and maintained below 10 KOhm by repositioning the electrodes as required. Using an inbuilt artefact rejection algorithm, the electrophysiology software automatically detected artefacts (e.g., from blinking), removed the corresponding responses and retested the sequence. Imaging, ERG, and optometric measurements were performed by certified technicians or ophthalmologists.

Patient diagnosis for DR or other eye issues was established once by experts at IMO

(M.G.R.), INDEREB (R.G.F.), APEC (H.Q.), and CSV (L.F.H.Z.) and after the finalization of our predictive diagnosis model, a second diagnosis (E.L.S. and V.C.S. from IMO) was established to compare our model performance with that of experts. The second diagnosis by experts was made under three conditions, the first with access to full eye examinations, the second with access to full eye tests and reference blood tests for diabetes, and the last with access to full eye tests, full blood tests, and anamnesis.

To test if ERG data from one sensor could improve of our model performance, a group of patients tested with the mydriasis-free RETeval ERG device was enrolled at APEC, INDEREB, and CSV from January to June 2022. **Table 3** contains an overview of these patient demographics and biometrics.

#### ERG data processing

For spectral analysis of both human and animal ERGs, signals were initially low-pass filtered at 1 kHz and high-pass filtered at 0.1 and 0.3 Hz for animals and humans, respectively. Recordings with large artifacts (which exceeded  $\pm 100 \,\mu\text{V}$ ) were removed. Recordings from the two eyes were independently analyzed. Raw ERG signals were normalized between -1 and +1.

ERGs show discontinuous activity, reason for which the wavelet (Morlet) transform was used to analyze them. Analysis was carried out with the MATLAB-based fieldtrip toolbox implementing the wavelet method <sup>100</sup>. The data were analyzed using custom-made MATLAB scripts (MATLAB R2018; MathWorks). Spontaneous human ERG signals were transformed within consecutive epochs of 10, 30, and 60 s. The corresponding number of windows was 30, 10, and 5, respectively. Time and spectral resolutions were 0.01 s and 0.05 Hz, respectively. The wavelet transform data were represented as scalograms or normalized

power spectra obtained by averaging the wavelet transform throughout the recording. The latter were subsequently grand averaged across all samples, animals, and patients for each condition. The standard error of the mean of the power spectra was calculated across animals/subjects.

In exploratory analysis, frequency points were initially considered between 0.1 or 0.3 Hz and 1 kHz (not shown) in animals and humans, respectively, and then refined in ranges where activity was detected (0.1-10 Hz and 0.3-40 Hz for animals and humans, respectively). Oscillatory potentials (OPs) were digitally isolated from the scotopic B-wave using a 100–500 Hz digital filter.

#### Modeling structure and development

Four common classification algorithms with built-in feature selection (Random Forest, deep learning and linear and radial support vector machines) were applied on human datasets using the open-source R package caret (version 6.0–73) for support vector machines and the H<sub>2</sub>O package for Random Forest and deep learning. Random Forest was applied on animal datasets using the H<sub>2</sub>O package. In all cases, the final datasets were randomly divided into training (80 % of observations) and testing (20 %) sets. Only validation dataset results were reported. Random Forest parameters were tuned as follows: ntrees = c(50, 70, 90, 100, 150, 200, 250, 300, 350, 400, 450, 500), max\_depth = c(9, 10, 11, 12, 13, 14, 15, 16, 17, 19, 20), min\_rows = c(1,2,3,4). We used the option in caret (precision, accuracy, sensitivity, and specificity) and H<sub>2</sub>O (ROC, ROC AUC, and confusion matrix) to return class probabilities for all classifiers.

#### Model performance analysis

Classes were balanced for all predictions. We therefore computed test performance metrics, including ROC curves <sup>101,102</sup>, AUC-ROC, accuracy, sensitivity, specificity, precision, negative predictive value (NPV) <sup>103</sup>, and confusion matrix <sup>104</sup>. The performance of our model was contrasted to the second patient diagnosis by experts using precision, recall, and F1-score.

#### Statistical analysis

Statistical analyses were performed using Matlab (Statistics and Machine Learning Toolbox). Data are reported as mean  $\pm$  s.e.m. or  $\pm$  95 % confidence interval. All data showed normal distribution and equal variance according to the D'Agostino-Pearson omnibus and Levene tests, respectively. Statistical significance was therefore determined either using unpaired ttest or, for multiple comparisons, using a mixed ANOVA and Bonferroni post-hoc. Human metrics were analyzed using the Welch ANOVA that accounts for variance heterogeneity. For explanatory statistical analyses, the following variables were considered: BMI, hip and waist circumferences, TG, HDL, blood pressure, fasting plasma glycemia, HbA1c, plasma insulin, LDL, VLDL, CT, creatinine, HOMA-I, atherogenic index, normalized power of 0.3 to 40 Hz oscillations with a frequency resolution of 0.05 Hz, and peak frequencies in the 0.3-2, 10-20, and 20-40 Hz bands. The variables with higher variance (CoV function in R) were selected, i.e., TG, CT, fasting glycemia, LDL, systolic blood pressure, age, body weight, hip circumference, HDL, VLDL, normalized power of 0.45, 0.4, 0.5, 0.55, 0.6 Hz activity, AUC of the 0.3-40 Hz band, normalized power of 0.65, 0.35, 0.7, 0.3, 0.75, 0.8 Hz activity, and AUC of the 20-40 Hz band, (ordered in descending order of variance) to perform PCA analysis. PCA panel and biplots of PC1, PC2, and PC3 were generated in R (princomp, biplot). Linear discriminant analyses were performed in R (library mass; code available in

the provided Github access). Cases with missing values were omitted.

To calculate the correlation between the spontaneous oscillation (SO) frequency components of major variance (detailed above), we computed a PCA score for each patient by combining linearly each variable coefficient (PC1) multiplied by its real value. Correlation analysis was done by calculating the Pearson's R coefficient (cor function in R).

The F1 score was calculated to rank our models for the shortest ERG window.

#### Code availability

We made use of several open-source libraries to conduct our experiments: *caret* (<u>https://topepo.github.io/caret/</u>) and H<sub>2</sub>O (<u>http://docs.h2o.ai/</u>), which provide implementations of individual model components. To facilitate improved reproducibility of our data analyses, the R code and documentation for the analysis are available online (<u>https://github.com/airetinopathydx/AIRetinopathyDx\_</u>).

#### References

- 1. Stem, M. S., Gardner, T. W. & Comer, G. M. 42: Management of Diabetic Retinopathy. in *Therapy for Diabetes Mellitus and Related Disorders* (2014).
- 2. Das, T. *et al.* Recently updated global diabetic retinopathy screening guidelines: commonalities, differences, and future possibilities. *Eye* **35**, 2685–2698 (2021).
- 3. Lawlor, M. S. *et al.* Cost of a group translation of the diabetes prevention program: Healthy living partnerships to prevent diabetes. *American Journal of Preventive Medicine* **44**, S381–S389 (2013).
- 4. Liu, J. *et al.* Health and Economic Impacts of the National Menu Calorie Labeling Law in the United States: A Microsimulation Study. *Circulation: Cardiovascular Quality and Outcomes* **13**, 309–319 (2020).
- 5. Kaur, G. *et al.* Diagnostic accuracy of tests for type 2 diabetes and prediabetes: A systematic review and meta-analysis. *PLOS ONE* **15**, e0242415 (2020).
- 6. Kanaya, A. M. & Venkat Narayan, K. M. Prevention of type 2 diabetes: Data from recent trials. *Primary Care Clinics in Office Practice* vol. 30 511–526 (2003).
- Buzsáki, G. & Draguhn, A. Neuronal olscillations in cortical networks. *Science* vol. 304 1926–1929 (2004).
- 8. Behar-Cohen, F. & Loewenstein, A. Current and Future Treatments for Diabetic Retinopathy. *Pharmaceutics* 14, 812 (2022).
- 9. Hyllienmark, L., Maltez, J., Dandenell, A., Ludvigsson, J. & Brismar, T. EEG abnormalities with and without relation to severe hypoglycaemia in adolescents with type 1 diabetes. *Diabetologia* **48**, 412–419 (2005).
- 10. Zeighami, Y. *et al.* Spontaneous neural activity changes after bariatric surgery: A resting-state fMRI study. *Neuroimage* **241**, (2021).
- 11. Biehl, S. C., Keil, J., Naumann, E. & Svaldi, J. ERP and oscillatory differences in overweight/obese and normal-weight adolescents in response to food stimuli. *Journal of Eating Disorders* **8**, 14 (2020).
- 12. Tammela, L. I. *et al.* Brain electrical activity during food presentation in obese binge-eating women. *Clinical Physiology and Functional Imaging* **30**, 135–140 (2010).
- 13. Imperatori, C. *et al.* Modification of EEG functional connectivity and EEG power spectra in overweight and obese patients with food addiction: An eLORETA study. *Brain Imaging and Behavior* **9**, 703–716 (2015).
- 14. Baver, S. B. *et al.* Leptin modulates the intrinsic excitability of AgRP/NPY neurons in the arcuate nucleus of the hypothalamus. *Journal of Neuroscience* **34**, 5486–5496 (2014).
- Pang, J. J., Gao, F. & Wu, S. M. Light-evoked excitatory and inhibitory synaptic inputs to ON and OFF α ganglion cells in the mouse retina. *Journal of Neuroscience* 23, 6063–6073 (2003).
- 16. Murphy, G. J. & Rieke, F. Network Variability Limits Stimulus-Evoked Spike Timing Precision in Retinal Ganglion Cells. *Neuron* **52**, 511–524 (2006).
- Margolis, D. J. & Detwiler, P. B. Different mechanisms generate maintained activity in ON and OFF retinal ganglion cells. *Journal of Neuroscience* 27, 5994–6005 (2007).
- 18. KUFFLER, S. W. Discharge patterns and functional organization of mammalian retina. *J Neurophysiol* **16**, 37–68 (1953).

- 19. Steinberg, R. H. Oscillatory activity in the optic tract of cat and light adaptation. *J Neurophysiol* **29**, 139–156 (1966).
- 20. Neuenschwander, S., Castelo-Branco, M. & Singer, W. Synchronous oscillations in the cat retina. Vision Research vol. 39 (1999).
- Petit-Jacques, J., Völgyi, B., Rudy, B. & Bloomfield, S. Spontaneous oscillatory activity of starburst amacrine cells in the mouse retina. *Journal of Neurophysiology* 94, 1770–1780 (2005).
- 22. Burrone, J. & Lagnado, L. *Electrical resonance and Ca2+ influx in the synaptic terminal of depolarizing bipolar cells from the goldfish retina. Journal of Physiology* (1997).
- 23. Ma, Y.-P. & Pan, Z.-H. Spontaneous regenerative activity in mammalian retinal bipolar cells: roles of multiple subtypes of voltage-dependent {Ca2}+ channels. *Visual Neuroscience* **20**, 131–139 (2003).
- 24. Solessio, E., Vigh, J., Cuenca, N., Rapp, K. & Lasater, E. M. Membrane properties of an unusual intrinsically oscillating, wide-field teleost retinal amacrine cell. *The Journal of Physiology* **544**, 831–847 (2002).
- Feigenspan, A., Gustincich, S., Bean, B. P. & Raviola, E. Spontaneous activity of solitary dopaminergic cells of the retina. *Journal of Neuroscience* 18, 6776–6789 (1998).
- 26. MATTHEWS, B. H. C. Sensory Mechanisms of the Retina. *Nature* **163**, 423–424 (1949).
- 27. Yee, C. W., Toychiev, A. H. & Sagdullaev, B. T. Network deficiency exacerbates impairment in a mouse model of retinal degeneration. *Frontiers in Systems Neuroscience* **6**, 8 (2012).
- 28. Trenholm, S. & Awatramani, G. B. Origins of spontaneous activity in the degenerating retina. *Frontiers in Cellular Neuroscience* **9**, 277 (2015).
- 29. Goo, Y. S., Park, D. J., Ahn, J. R. & Senok, S. S. Spontaneous Oscillatory Rhythms in the Degenerating Mouse Retina Modulate Retinal Ganglion Cell Responses to Electrical Stimulation. *Frontiers in Cellular Neuroscience* **9**, 512 (2016).
- 30. Euler, T. & Schubert, T. Multiple Independent Oscillatory Networks in the Degenerating Retina. *Frontiers in Cellular Neuroscience* **9**, 444 (2015).
- 31. Eggers, E. D. & Carreon, T. A. The effects of early diabetes on inner retinal neurons. *Visual Neuroscience* vol. 37 (2020).
- 32. Robson, A. G. *et al.* {ISCEV} guide to visual electrodiagnostic procedures. *Documenta Ophthalmologica. Advances in Ophthalmology* **136**, 1–26 (2018).
- 33. Kato, K., Kondo, M., Sugimoto, M., Ikesugi, K. & Matsubara, H. Effect of pupil size on flicker ERGs recorded with RETeval system: New mydriasis-free full-field ERG system. *Investigative Ophthalmology and Visual Science* **56**, 3684–3690 (2015).
- 34. Midena, E. *et al.* Early Microvascular and Oscillatory Potentials Changes in Human Diabetic Retina: Amacrine Cells and the Intraretinal Neurovascular Crosstalk. *Journal of Clinical Medicine* **10**, 4035 (2021).
- 35. Pinilla, I. *et al.* CHANGES IN TOTAL AND INNER RETINAL THICKNESSES IN TYPE 1 DIABETES WITH NO RETINOPATHY AFTER 8 YEARS OF FOLLOW-UP. *Retina* **40**, 1379–1386 (2020).
- 36. Pitale, P. M. & Gorbatyuk, M. S. Diabetic Retinopathy: From Animal Models to Cellular Signaling. *International Journal of Molecular Sciences* **23**, 1487 (2022).

- 37. Brigell, M. G., Chiang, B., Maa, A. Y. & Quentin Davis, C. Enhancing risk assessment in patients with diabetic retinopathy by combining measures of retinal function and structure. *Translational Vision Science and Technology* **9**, 1–9 (2020).
- 38. Goodfellow, I. J. *et al.* Generative Adversarial Networks. (2014) doi:10.48550/arxiv.1406.2661.
- 39. Festag, S., Denzler, J. & Spreckelsen, C. Generative adversarial networks for biomedical time series forecasting and imputation. *Journal of Biomedical Informatics* **129**, 104058 (2022).
- Hassan, W. M., Al-Dbass, A., Al-Ayadhi, L., Bhat, R. S. & El-Ansary, A. Discriminant analysis and binary logistic regression enable more accurate prediction of autism spectrum disorder than principal component analysis. *Scientific Reports* 12, 1–13 (2022).
- 41. Yuan, A. & Lee, A. Y. Artificial intelligence deployment in diabetic retinopathy: the last step of the translation continuum. *The Lancet Digital Health* vol. 4 e208–e209 (2022).
- 42. Mathenge, W. *et al.* Impact of artificial intelligence assessment of diabetic retinopathy on referral service uptake in a low resource setting: The RAIDERS randomized trial. *Ophthalmology Science* 100168 (2022) doi:10.1016/j.xops.2022.100168.
- 43. Bang, H. Development and Validation of a Patient Self-assessment Score for Diabetes Risk. *Annals of Internal Medicine* **151**, 775 (2009).
- 44. Pérez-Martínez, P. *et al.* Lifestyle recommendations for the prevention and management of metabolic syndrome: An international panel recommendation. *Nutrition Reviews* **75**, 307–326 (2017).
- 45. Pluta, W., Dudzińska, W. & Lubkowska, A. Metabolic Obesity in People with Normal Body Weight (MONW)—Review of Diagnostic Criteria. *International Journal of Environmental Research and Public Health* vol. 19 (2022).
- 46. Blüher, M. Metabolically healthy obesity. *Endocrine Reviews* vol. 41 405–420 (2020).
- 47. Simó, R., Stitt, A. W. & Gardner, T. W. Neurodegeneration in diabetic retinopathy: does it really matter? *Diabetologia* **61**, 1902–1912 (2018).
- 48. Keeling, E. *et al.* A High Fat "Western-style" Diet Induces AMD-like Features in Wildtype Mice. *Molecular Nutrition & Food Research* 2100823 (2022) doi:10.1002/mnfr.202100823.
- 49. Asare-Bediako, B. *et al.* Characterizing the Retinal Phenotype in the High-Fat Diet and Western Diet Mouse Models of Prediabetes. *Cells* **9**, 464 (2020).
- 50. Rajagopal, R. *et al.* Functional deficits precede structural lesions in mice with high-fat diet-induced diabetic retinopathy. *Diabetes* **65**, 1072–1084 (2016).
- 51. Lee, V. K. *et al.* BTBR ob/ob mouse model of type 2 diabetes exhibits early loss of retinal function and retinal inflammation followed by late vascular changes. *Diabetologia* **61**, 2422–2432 (2018).
- 52. King, A. J. The use of animal models in diabetes research. *British Journal of Pharmacology* **166**, 877–894 (2012).
- 53. Buettner, R., Schölmerich, J. & Bollheimer, L. C. High-fat diets: Modeling the metabolic disorders of human obesity in rodents. *Obesity* vol. 15 798–808 (2007).

- 54. Erkaymaz, O., Senyer Yapici, Í. & Uzun Arslan, R. Effects of obesity on timefrequency components of electroretinogram signal using continuous wavelet transform. *Biomedical Signal Processing and Control* **66**, 102398 (2021).
- 55. Yapici, İ. S., Erkaymaz, O. & Arslan, R. U. A hybrid intelligent classifier to estimate obesity levels based on ERG signals. *Physics Letters, Section A: General, Atomic and Solid State Physics* **399**, 127281 (2021).
- 56. Wachtmeister, L. Oscillatory potentials in the retina: What do they reveal. *Progress in Retinal and Eye Research* 17, 485–521 (1998).
- 57. Zelano, C. *et al.* Nasal respiration entrains human limbic oscillations and modulates cognitive function. *Journal of Neuroscience* **36**, 12448–12467 (2016).
- 58. Takahashi, K. & Atsumi, Y. Precise Measurement of Individual Rapid Eye Movements in REM Sleep of Humans. *Sleep* **20**, 743–752 (1997).
- 59. Mason, J. W. *et al.* Electrocardiographic reference ranges derived from 79,743 ambulatory subjects. *Journal of Electrocardiology* **40**, 228-234.e8 (2007).
- 60. Nair, G. *et al.* Effects of common anesthetics on eye movement and electroretinogram. *Documenta Ophthalmologica* **122**, 163–176 (2011).
- 61. Zenisek, D. & Matthews, G. Calcium action potentials in retinal bipolar neurons. *Visual Neuroscience* **15**, 69–75 (1998).
- 62. Vigh, J., Solessio, E., Morgans, C. W. & Lasater, E. M. Ionic mechanisms mediating oscillatory membrane potentials in wide-field retinal amacrine cells. *Journal of Neurophysiology* **90**, 431–443 (2003).
- 63. Trenholm, S. & Awatramani, G. B. *Myriad roles for gap junctions in retinal circuits*. *Webvision: The Organization of the Retina and Visual System* (University of Utah Health Sciences Center, 1995).
- 64. Hutcheon, B. & Yarom, Y. Resonance, oscillation and the intrinsic frequency preferences of neurons. *Trends in Neurosciences* vol. 23 216–222 (2000).
- Calbiague, V. M., Vielma, A. H., Cadiz, B., Paquet-Durand, F. & Schmachtenberg, O. Physiological assessment of high glucose neurotoxicity in mouse and rat retinal explants. *Journal of Comparative Neurology* 528, 989–1002 (2020).
- Nishimura, C. & Kuriyama, K. Alterations in the Retinal Dopaminergic Neuronal System in Rats with Streptozotocin-Induced Diabetes. *Journal of Neurochemistry* 45, 448–455 (1985).
- 67. Moore-Dotson, J. M. & Eggers, E. D. Reductions in calcium signaling limit inhibition to diabetic retinal rod bipolar cells. *Investigative Ophthalmology and Visual Science* **60**, 4063–4073 (2019).
- 68. Moore-Dotson, J. M. *et al.* Early retinal neuronal dysfunction in diabetic mice: Reduced light-evoked inhibition increases rod pathway signaling. *Investigative Ophthalmology and Visual Science* **57**, 1418–1430 (2016).
- Castilho, Á., Madsen, E., Ambrósio, A. F., Veruki, M. L. & Hartveit, E. Diabetic hyperglycemia reduces Ca<sup>2+</sup> permeability of extrasynaptic AMPA receptors in AII amacrine cells. *Journal of Neurophysiology* 114, 1545–1553 (2015).
- 70. Castilho, Á., Ambrósio, A. F., Hartveit, E. & Veruki, M. L. Disruption of a neural microcircuit in the rod pathway of the mammalian retina by diabetes mellitus. *Journal of Neuroscience* **35**, 5422–5433 (2015).
- 71. Yu, J. *et al.* Hyperactivity of ON-Type Retinal Ganglion Cells in Streptozotocin-Induced Diabetic Mice. *PLoS ONE* **8**, e76049 (2013).

- 72. Archana, R., Rajalakshmi, T. & Vijay Sai, P. Non-invasive technique to detect diabetic retinopathy based on Electrooculography signal using machine learning classifiers. *Proceedings of the Institution of Mechanical Engineers, Part H: Journal of Engineering in Medicine* 095441192210854 (2022) doi:10.1177/09544119221085422.
- 73. Wainio-Theberge, S., Wolff, A. & Northoff, G. Dynamic relationships between spontaneous and evoked electrophysiological activity. *Communications Biology* **4**, 1–17 (2021).
- Fan, F. L., Xiong, J., Li, M. & Wang, G. On Interpretability of Artificial Neural Networks: A Survey. *IEEE Transactions on Radiation and Plasma Medical Sciences* 5, 741–760 (2021).
- 75. Scuderi, G., Troiani, E. & Minnella, A. M. Gut Microbiome in Retina Health: The Crucial Role of the Gut-Retina Axis. *Frontiers in Microbiology* vol. 12 4246 (2022).
- Ting, D. S. W. *et al.* Development and validation of a deep learning system for diabetic retinopathy and related eye diseases using retinal images from multiethnic populations with diabetes. *JAMA - Journal of the American Medical Association* 318, 2211–2223 (2017).
- 77. FDA permits marketing of artificial intelligence-based device to detect certain diabetes-related eye problems | FDA. https://www.fda.gov/news-events/press-announcements/fda-permits-marketing-artificial-intelligence-based-device-detect-certain-diabetes-related-eye.
- 78. Burke, H. B. The power of prediction. *Cancer* vol. 113 890–892 (2008).
- 79. Gravina, R., Alinia, P., Ghasemzadeh, H. & Fortino, G. Multi-sensor fusion in body sensor networks: State-of-the-art and research challenges. *Information Fusion* **35**, 1339–1351 (2017).
- 80. Robson, A. G. *et al.* ISCEV Standard for full-field clinical electroretinography (2022 update). *Documenta Ophthalmologica* (2022) doi:10.1007/s10633-022-09872-0.
- 81. White, N. H. *et al.* The Effect of Interventions to Prevent Type 2 Diabetes on the Development of Diabetic Retinopathy: The DPP/DPPOS Experience. *Diabetes Care* (2022) doi:10.2337/dc21-2417.
- 82. Harrison, W. W. *et al.* Multifocal electroretinograms predict onset of diabetic retinopathy in adult patients with diabetes. *Investigative Ophthalmology and Visual Science* **52**, 772–777 (2011).
- 83. BEST (Biomarkers, EndpointS, and other Tools) Resource [Internet] PubMed. https://pubmed.ncbi.nlm.nih.gov/27010052/.
- 84. Burke, H. B. Predicting Clinical Outcomes Using Molecular Biomarkers. *Biomarkers in Cancer* **8**, BIC.S33380 (2016).
- 85. Sachdeva, M. M. Retinal Neurodegeneration in Diabetes: an Emerging Concept in Diabetic Retinopathy. *Current Diabetes Reports* vol. 21 (2021).
- 86. The increasing potential of electroretinography in ophthalmic diagnostics. https://europe.ophthalmologytimes.com/view/the-increasing-potential-ofelectroretinography-in-ophthalmic-diagnostics.
- 87. Prior, L. *et al.* Patients views on a new surveillance pathway involving allied nonmedical staff for people with treated diabetic macular oedema and proliferative diabetic retinopathy. *Eye* 1–5 (2022) doi:10.1038/s41433-022-02050-1.
- 88. Miranda-Anaya, M., Pérez-Mendoza, M., Juárez-Tapia, C. R. & Carmona-Castro, A. The volcano mouse {Neotomodon} alstoni of central {Mexico}, a biological model

in the study of breeding, obesity and circadian rhythms. *General and Comparative Endocrinology* **273**, 61–66 (2019).

- 89. Like, A. A. & Rossini, A. A. Streptozotocin-induced pancreatic insulitis: new model of diabetes mellitus. *Science* **193**, 415–417 (1976).
- 90. King, A. J. F. The use of animal models in diabetes research. *British Journal of Pharmacology* **166**, 877–894 (2012).
- 91. Han, B. G. *et al.* Markers of glycemic control in the mouse: comparisons of 6-h- and overnight-fasted blood glucoses to {Hb} {A1c}. *American Journal of Physiology. Endocrinology and Metabolism* **295**, E981--986 (2008).
- 92. Improved contact lens electrode for corneal {ERG} recordings in mice {\textbar} {SpringerLink}.
- 93. Armington, J. C. Electroretinogram, {Electroretinography}. in Sensory {System} {I}: {Vision} and {Visual} {Systems} (ed. Held, R.) 21–22 (Birkhäuser, 1988). doi:10.1007/978-1-4899-6647-6\_11.
- 94. Matthews, D. R. *et al.* Homeostasis model assessment: insulin resistance and betacell function from fasting plasma glucose and insulin concentrations in man. *Diabetologia* **28**, 412–419 (1985).
- 95. Global report on diabetes.
- 96. Alberti, K. G. & Zimmet, P. Z. Definition, diagnosis and classification of diabetes mellitus and its complications. {Part} 1: diagnosis and classification of diabetes mellitus provisional report of a {WHO} consultation. *Diabetic Medicine: A Journal of the British Diabetic Association* **15**, 539–553 (1998).
- 97. Consensus statements.
- 98. Nwagha, U. I., Ikekpeazu, E. J., Ejezie, F. E., Neboh, E. E. & Maduka, I. C. Atherogenic index of plasma as useful predictor of cardiovascular risk among postmenopausal women in {Enugu}, {Nigeria}. *African Health Sciences* **10**, 248–252 (2010).
- 99. Zeng, Y. *et al.* Screening for diabetic retinopathy in diabetic patients with a mydriasis-free, full-field flicker electroretinogram recording device. *Documenta Ophthalmologica. Advances in Ophthalmology* **140**, 211–220 (2020).
- Oostenveld, R., Fries, P., Maris, E. & Schoffelen, J.-M. {FieldTrip}: {Open} source software for advanced analysis of {MEG}, {EEG}, and invasive electrophysiological data. *Computational Intelligence and Neuroscience* 2011, 156869 (2011).
- 101. DeVries, P. M. R., Viégas, F., Wattenberg, M. & Meade, B. J. Deep learning of aftershock patterns following large earthquakes. *Nature* **560**, 632–634 (2018).
- 102. Provost, F., Fawcett, T. & Kohavi, R. The {Case} {Against} {Accuracy} {Estimation} for {Comparing} {Induction} {Algorithms}. in *In {Proceedings} of the {Fifteenth} {International} {Conference} on {Machine} {Learning} 445–453* (Morgan Kaufmann, 1998).
- 103. Neural network-based approaches for biomedical relation classification: {A} review
  {ScienceDirect}.
- 104. Dankers, F. J. W. M., Traverso, A., Wee, L. & van Kuijk, S. M. J. Prediction {Modeling} {Methodology}. in *Fundamentals of {Clinical} {Data} {Science}* (eds. Kubben, P., Dumontier, M. & Dekker, A.) (Springer, 2019).

#### Acknowledgments

We thank all volunteers under the care of the AEPC, IMO, and INDEREB. We thank Edgar Morales, Oliver Becerra, Javier Ledezma, Marco Arieli Herrera Valdez, Adriana Petriz, Susan Pérez Salazar, Jesús Chávez Baldera, and Dulce María Soria Lara for their participation in data acquisition and/or processing. We thank Edith Espino, Marina Ramírez Romero, Martín García, and Alejandra Castilla for their technical assistance, and Jessica Norris for critically editing the manuscript. We thank the many physicians, technicians, and administrators of the IMO, APEC, INDEREB, and ENES León who worked on the data collection. We thank Roger D. Traub M.D. for helpful discussions, as well as Pavel Rueda Orozco Ph.D. and Hugo Merchant Ph.D. for insightful comments. R.N.I. is a Doctoral student from the Programa de Posgrado en Ciencias, Universidad Nacional Autónoma de México (UNAM) and received fellowships from the National Council of Science and Technology of Mexico (CONACYT; #620199). This study was supported by the UNAM grant IN209317 (ST), IN205420 (ST), CONACYT 299625 (ST), CONACYT CF-2019-1759 (ST), and Shedid grant (R. Miledí and A. Martínez Torres).

S.C.T. conceived and designed the experiments, analyzed the data, wrote the first draft of the manuscript, and contributed to the writing of the manuscript. J.-M.B., D.M., E.B., and G.-M.C. developed the random forest-based models and performed model analysis. R.N.I. performed *in vivo* ERGs in experimental models, analyzed the data, and prepared figures. J.-M-B. developed GAN model and created GAN data. M.M.A. provided spontaneously obese mice and metabolic data. P.R.O. performed PCA, cluster, and dendogram analysis under the supervision of L.G.-M. H.Q.-M-, L.F.H-Z., J.A.H.C., R.G.-F., E.L.-S., M.G.R., V.C.L. provided human data and supervised ophthalmological diagnosis of patients. S.C.T., F.P.-O.,

G.R.-P. agree with manuscript results and conclusions, jointly developed some arguments for the paper, and made critical revisions. All authors approved the final manuscript.

The authors declare that they have no competing interests.

All experimental data needed to evaluate the conclusions in the paper are present in the paper. The clinical data used for the training and test sets were collected at the IMO, APEC, INDEREB, and ENES León and transferred to a secure data center with strict access controls in de-identified format. Clinical data were used with individual permissions. They are not publicly available, and restrictions apply for their use. The de-identified dataset (or a test subset) can be provided by Stéphanie C. Thébault pending scientific review and a completed material transfer agreement. Requests for the clinical data should be submitted to: <a href="mailto:sthebault@comunidad.unam.mx">sthebault@comunidad.unam.mx</a>.

#### **Figure legends**

Figure 1: Spontaneous ERG-based random forest model discriminates control and disease cases in rodent models of obesity and type 1 diabetes and predicts disease evolution in the type 1 diabetes model. Illustrative spontaneous ERGs and wavelet analysis in A, control versus high-fat diet-fed mice (n = 75 and n = 75, respectively), C, lean versus spontaneously obese *Neotomodon alstoni* mice (n = 20 and n = 20, respectively), and E, control and streptozotocin-treated rats (n = 40 and n = 40, respectively) in the 0.1-10 Hz range, under photopic conditions. ERG signals were normalized. Graphs show the average scalogram power  $\pm$  s.e.m. throughout 1-minute recordings. The square brackets with the Roman numerals indicate the consistent peaks observed in each control condition. In the high-fat diet-induced obesity model, 0.1-0.8 (Low, L), 1-1.8 (M1, mid-low), and 2-4 Hz (M2, mid) bands were considered; 0.1-0.6 (I), 0.6-1 (II), and 1-1.7 (III) Hz bands in the spontaneous model of obesity, and 0.2-0.6 (I) and 0.6-2.5 (II) Hz bands in the type 1 diabetes model. AUC (0.1-10 Hz) and peak frequency analysis (P values were determined by unpaired Student's t-test) of wavelet graphs in **B**, the high-fat diet-induced obesity **D**, the spontaneous model of obesity, and F, the type 1 diabetes models. Graphs show mean  $\pm$  confidence interval. ROC curves and confusion matrix with performance measures for binary predictions (control vs. experimental cases) using the 0.1-10 Hz power spectra of G, control and highfat diet-fed mice (n = 15 and n = 15, respectively), H, lean and spontaneously obese Neotomodon alstoni mice (n = 4 and n = 4, respectively), and I, control and streptozotocintreated rats (n = 8 and n = 8, respectively). S, sensitivity; Sp., Specificity; NPV, negative predictive value; Ac., accuracy. J, Confusion matrix of multiclass prediction for the machine learning algorithm that discriminates between control and diseased rats at week 4, 6, 8, or 12 post-streptozotocin injection. Each column represents the instances in a predicted class, and

rows represent the instances in an actual class. We used controls at week 4, 6, 8, and 12 (n = 8 at each time point).

Figure 2: Spontaneous ERG oscillations help predict the disease in patients with early risk factors of DR. A, Illustrative signals from spontaneous ERGs of control subjects and patients with overweight (OW), obesity, metabolic syndrome (MetS), and diabetes but no DR (DM no DR) under photopic conditions. Signals are normalized. Wavelet analysis for **B**, each group separately and C, combining all unhealthy conditions into one disease group. Control group is defined in Methods. Graphs show the average scalogram power  $\pm$  s.e.m. throughout 20-second recordings between 0.3-40 Hz (host graph) and 0.3-2 (I), 10-20 (II), and 20-40 (III) Hz (inset graphs), where consistent peaks were observed in the control group. **D**, AUC and peak frequency analysis in the I, II, and III bands of wavelet power spectra from control and disease groups (P values were determined by unpaired Student's t-test). Graphs show mean  $\pm$  confidence interval. **E**, ROC curve and **F**, confusion matrix with performance measures for the Random Forest model discriminating control from disease cases in control and disease groups, using the 0.3-40 Hz power spectra. G, ROC curve and H, confusion matrix with performance measures corresponding to the validation of our predictive model thanks to an external and independent validation dataset. I, Classification performance of our model versus experts, when experts have access to full eye examination alone or combined with reference tests to diagnose diabetes (fasting glycemia, creatinine, triglycerides, and cholesterol) or full laboratory (glycemia, HbA1c levels, HOMA-I, and lipid profile) and nonlaboratory (blood pressure, weight and waist circumference, and BMI) tests. Precision, recall, and F1-score are reported. I, ROC curve and J, confusion matrix with performance measures corresponding to the performance of our predictive model using ERG-derived power spectra

recorded by only one ERG device. In **A-F** and **I**, control (metabolically healthy, n = 80) and disease (OW, n = 40; obese, n = 14; MetS, n = 66; and diabetes with no DR, n = 68) groups. In **G**, **H**, control (metabolically healthy, n = 11) and disease (OW, n = 5; obese, n = 3; MetS, n = 3; and diabetes with no DR, n = 3) groups. In **J**, **K**, control (metabolically healthy, n =37) and disease (OW, n = 18; obese, n = 8; MetS, n = 42; and diabetes with no DR, n = 33) groups. S, sensitivity; Sp., Specificity; NPV, negative predictive value; Ac., accuracy.

Figure 3: Multiclass prediction of people at risk of developing type 2 diabetes and DR using the spontaneous ERG-based random forest model: explanatory data analysis. A, Multiclass confusion matrix for the prediction of control, overweight (OW), obesity, MetS, and diabetes with no DR (DM no DR) once the first binary classification (control vs. disease) was done using the Random Forest model. Data from Table 1 and 2 were pooled and aleatory divided in a training (80 %) and test (20 %) set. **B**, Separation of control, overweight, obese, metabolic syndrome, and diabetes (DM) without DR cases using PCA and C, corresponding biplot of PC1 and PC2, using 13 variables: normalized power of the slowest oscillations (0.3-0.8 Hz), the AUC of the 0.3-40 Hz band (0.3-40 Hz AUC), and the peak frequency of the 20-40 Hz band (20-40 Hz PF). D, Separation of control, overweight, obese, metabolic syndrome, and diabetes (DM) without DR cases using PCA and E, corresponding biplot of PC1 and PC2, using 20 variables: normalized power of 0.4-0.7 Hz activities, TG, CT, fasting blood glucose (BG), LDL, systolic blood pressure (SBP), VLDL, age, body weight (BWT), hip and waist circumferences (HC and WC, respectively), HDL, insulinemia (Ins.), and diastolic blood pressure (DBP). Table (right) shows variable contributions to the first and second principal components. F, Separation of control, overweight, obese, metabolic syndrome, and diabetes (DM) without DR cases using PCA, using 8 variables: normalized power of 0.7 to

0.8 Hz activities, the peak frequency in the 20-40 Hz band, the AUC of the 0.3-40 Hz band, age, BG, CT. **G**, Confusion matrix corresponding to the linear discriminant analysis to separate control, diabetes risk factors (overweight, obesity, and MetS, n = 168), and diabetes (DM) without DR groups. From **B** to **G**, control (n = 62), overweight (OW, n = 66), obese (n = 17), metabolic syndrome (MetS, n = 85), and diabetes (DM) without DR (n = 41).

Supplementary Figure 1: Metabolic follow-up of animal models and oscillatory potential analysis. A, Follow-up of body weight and blood glucose levels, glucose tolerance test, and insulin tolerance test in control-diet (n = 75) and high-fat diet-fed (n = 75) mice for 12 weeks. **B**, Body weight, blood glucose levels, glucose tolerance test, and insulin tolerance test in lean (n = 20) and spontaneously obese (n = 20) Neotomodon alstoni mice. In both models, mice have higher glycemia than control mice at every time point (P < 0.05), suggesting reduced insulin sensitivity. This was confirmed by insulin tolerance tests, which showed a lower fall in blood glucose in response to insulin in obese mice as compared with control mice. C, Blood glucose level follow-up in rats after 4, 6, 8, or 12 weeks of streptozotocin (n = 40) and vehicle (n = 40) treatment. Values, mean  $\pm$  s.d. \* indicates P values < 0.05 determined by a two-sample Student's *t*-test in **B** (body weight and glycemia) and by a mixed ANOVA followed by Bonferroni test everywhere else. D, Illustrative ERG (top) and oscillatory potentials (OP, bottom) in control-diet and high-fat diet-fed mice for 12 weeks, measured in response to a light flash of 7.72 (cd.s)/m<sup>2</sup> (arrow) under dark-adapted conditions. E, Temporal monitoring (0 to 12 weeks, as indicated) of the average amplitude and implicit time of OP1, OP2, OP3, and OP4 in control (n = 10) and high-fat diet fed mice (n = 12) under dark-adapted conditions at increasing light intensities (0.02, 0.24, 2.45, and 7.72 (cd.s)/m<sup>2</sup>. CD, control diet. HFD, high-fat diet. n.s., not significant (P > 0.05).

Supplementary Figure 2: Oscillatory potentials and DR score analysis in people with early risk factors for DR, and correlation with slow frequency spontaneous ERG oscillation power. A, Illustrative recording from the ISCEV DA 3 ERG protocol in a control case. N1-P1 and P1-N2 B, amplitude and C, peak time ratio in control (metabolically healthy, n = 8) and disease (OW, n = 5; obese, n = 2; MetS, n = 6; and diabetes with no DR, n = 61) groups. D, DR score in control (metabolically healthy, n = 74) and disease (OW, n = 22; obese, n = 18; MetS, n = 88; and diabetes with no DR, n = 109) groups. B-D, Graphs show mean  $\pm$  confidence interval. Correlation analysis between the E, N1-P1 or P1-N2 amplitude ratio, F, N1-P1 or P1-N2 peak time ratio, and G, DR score with the spontaneous oscillation (SO) PCA score (detailed in Methods) in the groups of interest.

Supplementary Figure 3: Random Forest model performs better than support vector machine algorithms and when it primarily uses spontaneous photopic ERG of 60-s duration in humans. A, ROC curves for both linear and radial svm algorithms. B, Performance parameters for the random forest model using power spectra from photopic or mesopic ERGs of 10, 30 or 60 s. C, ROC curves for the random forest model using power spectra from photopic, mesopic or combined photopic and mesopic ERGs of 60 s. D, Corresponding performance parameters. All data correspond to binary classification between control and disease cases. Controls are constituted by metabolically healthy subjects (n = 62) and the disease group by patients with overweight (n = 41), obesity (n = 16), metabolic syndrome (n = 55), and diabetes with no DR (n = 63).

Table 1: Characteristics of the patients studied for model training and test using ERG spectral components from three different recording devices. Values, mean  $\pm$  s.d. *P* values were determined by the Welch ANOVA test. Values that do not share a letter (a, b, or c) are statistically different. Years, y.

Table 2: Characteristics of the patients studied for external model validation. Values, mean  $\pm$  s.d. *P* values were determined by the Welch ANOVA test. Values that do not share a letter (a, b, or c) are statistically different. Years, y.

Table 3: Characteristics of the patients studied for external model validation using ERG spectral components from one recording device. Values, mean  $\pm$  s.d. *P* values were determined by the Welch ANOVA test. Values that do not share a letter (a, b, or c) are statistically different. Years, y.















PC	C1 (26%)	PC2 (46%)
TG	0.092	0.304
CT	0.064	0.118
BG	0.065	0.241
LDL	0.054	0.071
SBP	0.028	0.321
VLDL	0.107	0.289
Age	-0.063	0.200
BWT	0.117	0.332
HC	0.092	0.355
WC	0.108	0.333
HDL	-0.047	-0.228
Ins.	0.069	0.199
DBP	0.106	0.305
0.5	-0.366	0.091
0.45	-0.358	0.058
0.55	-0.390	0.124
0.6	-0.385	0.132
0.65	-0.370	0.125
0.7	-0.341	0.106
0.4	-0.311	0.034



F •Control •Obese •DM no DR •OW •MetS



G			
	Control	Diabetes	Diabetes RF
Control	47	3	12
Diabetes	6	39	21
Diabetes RF	19	12	109

### Supplementary FIG. 1





kappa

-0.016

0.274

## Supplementary FIG. 3



С

Model	AUC	Accuracy	Precision	Error	Sensitivity	Specificity			
Disease (N = 35) vs. control (N = 12)									
Deep learning PHOTOPIC 10 s	0.670	0.55	0.82	0.45	0.21	0.94			
Deep learning PHOTOPIC 30 s	0.613	0.5	1	0.5	0.07	1			
Deep learning PHOTOPIC 60 s	0.809	0.77	0.83	0.23	0.71	0.83			
Deep learning MESOTOPIC 10 s	0.593	0.58	0.73	0.42	0.24	0.91			
Deep learning MESOTOPIC 30 s	0.606	0.6	1	0.4	0.19	1			
Deep learning MESOTOPIC 60 s	0.562	0.6	1	0.4	0.19	1			





Ε

Model	AUC	Accuracy	Precision	Error	Sensitivity	Specificity		
Disease (N = 28/7) vs. control (N = 24/6)								
Deep learning PHOTOPIC 60 s	0.809	0.77	0.83	0.23	0.71	0.83		
Deep learning MESOTOPIC 60 s	0.562	0.65	0.83	0.35	0.38	0.92		
Deep learning PHOTO/MESO 60 s	0.595	0.54	1	0.46	0.15	1		

## TABLE 1

	CONTROL	OVERWEIGHT	OBESITY	MetS	DIABETES Without RD	Ρ
n	80	40	14	66	68	
Age (y)	$41.49 \pm 1.88^{\text{d}}$	$42.70 \pm 2.08^{cd}$	39.57 ± 2.93 <sup>cd</sup>	60.45 ± 1.79 <sup>bc</sup>	65.86 ± 1.49ª	<0.0001
DM1	-	-	-	-	7	
DM2	-	-	-	-	61	
Duration of diabetes (y)	-	-	-	-	9.45 ± 1.04°	0.0077
Body weight (Kg)	58.38 ± 1.23°	72.47 ± 1.40 <sup>ab</sup>	79.12 ± 3.97ª	69.65 ± 4.41ª	76.48 ± 1.86ª	<0.0001
Waist circunference (cm)	83.13 ± 1.30°	96.05 ± 1.38 <sup>ab</sup>	104.14 ± 4.26ª	88 ± 3.77ª	102.91 ± 1.53ª	<0.0001
Abdominal circunference (cm)	82.28 ± 1.25 <sup>c</sup>	93.80 ± 1.38 <sup>b</sup>	99.46 ± 3.94ªb	96.36 ± 3.96 <sup>ab</sup>	103.17 ± 1.18ª	<0.0001
BMI (Kg/m²)	$22.39\pm0.31^{\text{d}}$	26.85 ± 0.25 <sup>bc</sup>	29.93 ± 1.82 <sup>abc</sup>	27.12 ± 1.08 <sup>ab</sup>	30.14 ± 0.63ª	<0.0001
Glycemia (mg/dl)	86.77 ± 0.79 <sup>b</sup>	87.75 ± 1.21 <sup>b</sup>	91.43 ± 3.09 <sup>b</sup>	144.82 ± 20.1 <sup>b</sup>	147.48 ± 6.39ª	<0.0001
HbA1c (%)	$5.30\pm0.03^{d}$	$5.47 \pm 0.05^{d}$	$5.39 \pm 0.10^{d}$	$8.41\pm0.69^{\text{d}}$	7.67 ± 0.24 <sup>c</sup>	<0.0001
Insulinemia (μUI/ml)	7.36 ± 0.44 <sup>c</sup>	9.54 ± 0.91°	11.15 ± 1.74 <sup>abc</sup>	$9.84 \pm 1.70^{bc}$	$21.50 \pm 2.89^{ab}$	<0.0001
HOMA-I	$1.67 \pm 0.10^{b}$	$2.06 \pm 0.21^{b}$	$2.64 \pm 0.49^{ab}$	$3.44 \pm 0.70^{b}$	8.66 ± 1.31ª	<0.0001
TG	100.37 ± 4.79ª	104. 70 ± 6.99 <sup>c</sup>	104.94 ± 12.56 <sup>bc</sup>	132.55 ± 21.4ª	$180.6 \pm 11.05^{ab}$	<0.0001
CT (mg/dl)	$179.55 \pm 3.61^{b}$	185.12 ± 4.37 <sup>ab</sup>	192.19 ± 8.05 <sup>ab</sup>	$190.17 \pm 12.66^{a}$	$183.36\pm4.39^{\text{b}}$	0.0015
HDL (mg/dl)	56.92 ± 2.02ª	54.64 ± 1.75 <sup>ab</sup>	$56.90 \pm 3.36^{abc}$	48.84 ± 5.21°	46.13 ± 1.45 <sup>c</sup>	<0.0001
LDL (mg/dl)	105 ± 2.72 <sup>b</sup>	117.27 ± 3.73 <sup>ab</sup>	118.25 ± 7.97 <sup>ab</sup>	113.06 ± 11.02ª	$108.55 \pm 3.84^{ab}$	0.00889
VLDL (mg/dl)	$17.23 \pm 1.08^{d}$	$17.62 \pm 1.98^{d}$	$17.05 \pm 1.94^{cd}$	26.51 ± 4.27 <sup>ab</sup>	30.34 ± 2.37 <sup>bc</sup>	<0.0001
Creatinine (mg/dl)	$0.79 \pm 0.02^{b}$	$0.77 \pm 0.02^{ab}$	$0.75 \pm 0.04^{\text{ab}}$	$1.09\pm0.18^{\rm ab}$	$0.79 \pm 0.05^{b}$	0.00965
Systolic blood pressure (mmHg)	114.73 ± 2.23°	114.79 ± 2.82°	124.14 ± 7.75 <sup>bc</sup>	125.64 ± 4.72 <sup>ab</sup>	131.33 ± 2.62 <sup>b</sup>	<0.0001
Diastolic blood pressure (mmHg)	$71.64 \pm 1.16^{d}$	$74.05 \pm 1.65^{cd}$	73.14 ± 3.21 <sup>bcd</sup>	85.82 ± 4.31ª	82.61 ± 1.50ªb	<0.0001

## TABLE 2

	CONTROL	OVERWEIGHT	OBESITY	MetS	DIABETES without RD	Р
n	19	6	4	5	4	
Age (y)	30.50 ± 3.12°	38.80 ± 7.62°	33 ± 2.12 <sup>bc</sup>	31.5 ± 2.47°	67.33 ± 7.31ª	< 0.0001
DM1	-	-	-	-	-	-
DM2	-	-	-	-	3	
Duration of diabetes (y)	-	-	-	-	2.55 <b>±</b> 1.73ª	<0.0001
Body weight (Kg) Waist	53.57 <b>±</b> 2.32 <sup>♭</sup>	70.98 ± 2.95ª	67.55 ± 1.52 <sup>ab</sup>	74.60 ± 3.25 <sup>a</sup>	82.20 ± 1.98 <sup>a</sup>	<0.0001
circunference (cm) Abdominal	86.13 ± 2.82 <sup>b</sup>	102.52 ± 2.23ª	89 ± 0.53ªb	101.75 ± 0.53 <sup>ab</sup>	95.80 ± 1.56 <sup>ab</sup>	<0.0001
circunference (cm)	73.56 <b>±</b> 3.75℃	88.94 ± 5.01 <sup>bc</sup>	97 ± 1.41 <sup>abc</sup>	93 ± 1.77 <sup>abc</sup>	105.90 ± 0.64 <sup>ab</sup>	<0.0001
BMI (Kg/m²)	20.38 ± 0.67°	26.40 ± 0.45 <sup>ab</sup>	26.66 ± 0.73 <sup>bc</sup>	$28.08 \pm 0.48^{\mathrm{ab}}$	27.97 ± 1.02 <sup>ab</sup>	<0.0001
Glycemia (mg/dl)	84.89 ± 2 <sup>b</sup>	83.60 ± 1.71 <sup>b</sup>	87 ± 5.66 <sup>b</sup>	89 ± 2.83 <sup>b</sup>	118 ± 4.50 <sup>b</sup>	<0.0001
HbA1c (%)	5.32 ±0.03ª	5.62 ± 0.12 <sup>a</sup>	5.20 ± 0.10 <sup>a</sup>	5.35 ± 0.32ª	$7.03 \pm 0.38^{b}$	<0.0001
Insulinemia (µUI/ml)	7.49 <b>±</b> 1.25 <sup>b</sup>	$10.14 \pm 1.84^{b}$	15.25 ± 6.26 <sup>b</sup>	$14.15 \pm 1.38^{b}$	18.30 ± 3.60 <sup>a</sup>	0.00245
HOMA-I	$1.61 \pm 0.29^{b}$	$1.84 \pm 0.26^{b}$	3.45 ± 1.59 <sup>b</sup>	$2.85 \pm 0.18^{b}$	5.40 ± 0.98 <sup>a</sup>	0.0483
TG	76.75 ± 8.74 <sup>b</sup>	93.84 ± 13.38 <sup>b</sup>	110.55 ± 1.47 <sup>b</sup>	290 ± 46.64ª	142.85 ± 17.1 <sup>b</sup>	<0.0001
CT (mg/dl)	177.21 <b>±</b> 4.26	173.13 ± 0.25	183.27 ± 0.07	205.56 ± 1.02	184.13 ± 0.32	0.5103
HDL (mg/dl)	55.69 ± 7.30 <sup>b</sup>	55.74 ± 3.20 <sup>b</sup>	45.49 ± 3.53ª	39.16 ± 6.45ª	44.36 ± 3.28 <sup>a</sup>	0.0345
LDL (mg/dl)	98.35 ± 7.66 <sup>b</sup>	112.26 ± 0.18 <sup>ab</sup>	111.94 ± 0.09 <sup>ab</sup>	123.27 ± 0.61ªb	102.6 ± 0.32 <sup>ab</sup>	0.04646
VLDL (mg/dl)	9.36 ± 1.57 <sup>bc</sup>	12.94 ± 2.52 <sup>bc</sup>	14.75 ± 0.51 <sup>bc</sup>	45.10 ± 7.77ª	33.33 ± 2.85 <sup>b</sup>	<0.0001
Creatinine (mg/dl)	$0.78 \pm 0.04^{\text{b}}$	$0.77 \pm 0.04^{b}$	$0.78 \pm 0.02^{b}$	$0.80 \pm 0.04^{ab}$	0.91 ± 0.08 <sup>ab</sup>	0.00418
Systolic blood pressure (mmHg)	106.89 ± 3.71	109 ± 6.41	110.50 ± 8.13	138.50 ± 3.18	131 ± 1.41	0.0643
pressure (mmHg)	73.33 <b>±</b> 3.43	76.80 ± 4.47	66 <b>±</b> 6.36	78 ± 3.54	86 ± 9.19	0.059

### TABLE 3

					DIABETES	
	CONTROL	OVERWEIGHT	OBESE	MetS	Without RD	р
n	37	18	8	42	33	
Age (y)	34.92 <b>±</b> 2.85℃	38.83 ± 3.59 <sup>bc</sup>	38.25 ± 3.12 <sup>bc</sup>	46.12 ± 2.02 <sup>b</sup>	60.58 ± 2.05ª	< 0.0001
DM1	-	-	-	-	1	
DM2	-	-	-	-	32	
Duration of						
diabetes (y)	-	-	-	-	6.15 ± 1.12 <sup>b</sup>	0.0003
Body weight (Kg)	55.03 <b>±</b> 1.45°	71.41 ± 1.80 <sup>ab</sup>	75.34 ± 6.62 <sup>ab</sup>	78.20 ± 2.03ª	72.19 ± 1.99 <sup>ab</sup>	<0.0001
Waist circunference (cm) Abdominal	83.24 ± 3.05 <sup>b</sup>	96.95 ± 1.80ª	97.81 ± 4.65ª	102.63 ± 1.30ª	102.25 ± 1.74ª	<0.0001
circunference						
(cm)	73.82 ± 2.94 <sup>b</sup>	86.69 ± 5.64 <sup>ab</sup>	96.28 ± 5.85°	98.92 ± 1.78 <sup>a</sup>	101.69 ± 1.41 <sup>a</sup>	<0.0001
BMI (kg/m²)	21.15 ± 0.35°	$26.78 \pm 0.34^{ab}$	26.48 ± 1.91 <sup>ab</sup>	29.03 ± 0.65ª	29.71 ± 0.73ª	<0.0001
Glycemia (mg/dl)	82.84 ± 1.08°	85.18 ± 1.14°	90.50 ± 1.85°	89.93 ± 1.87°	132.42 ± 8.56 <sup>b</sup>	< 0.0001
HbA1c (%)	5.32 ± 0.04 <sup>d</sup>	5.48 ± 0.06 <sup>d</sup>	5.31 ± 0.08 <sup>d</sup>	5.62 ± 0.06 <sup>d</sup>	7.42 ± 0.30°	< 0.0001
Insulinemia						
(µl/ml)	7.82 ± 0.64 <sup>b</sup>	9.99 ± 1.35 <sup>b</sup>	10.89 ± 1.92 <sup>b</sup>	13.65 ± 1.29 <sup>ab</sup>	18.22 ± 4.15 <sup>ab</sup>	0.001
HOMA-I	1.61 ± 0.12 <sup>b</sup>	2.02 ± 0.27 <sup>b</sup>	2.46 ± 0.47 <sup>b</sup>	3.11 ± 0.33 <sup>b</sup>	5.69 ± 1.12 <sup>ab</sup>	< 0.0001
TG	107.41 ± 8.26 <sup>b</sup>	96.89 ± 8.67 <sup>b</sup>	132.63 ± 18.93 <sup>b</sup>	210.66 ± 18.48 <sup>b</sup>	188.76 ± 13.18 <sup>b</sup>	<0.0001
CT (mg/dl)	176 20 <b>+</b> 4 53 <sup>b</sup>	171 70 + 7 88 <sup>b</sup>	186 57 + 11 43ab	201 50 + 5 56ab	182 93 + 6 09 <sup>b</sup>	0.012
HDI (mg/dl)	54.29 + 2.78ª	$53.09 + 2.95^{\circ}$	54 96 + 5,99ª	$44.62 + 1.41^{\circ}$	46.83 + 1.58ª	0.028
	54.25 2 2.70	<u> 55.05 ± 2.55</u>	54.50 ± 5.55	44.02 <u>-</u> 1.41	40.05 ± 1.50	0.020
LDL (mg/dl)	106.62 ± 3.70ª	113.24 <b>±</b> 6.35ª	112.39 ± 11.96ª	124.43 <b>±</b> 4.70ª	106.35 ± 5.33ª	0.128
VLDL (mg/dl)	14.84 ± 1.42 <sup>c</sup>	13.11 ± 1.50°	19.15 ± 3.22 <sup>bc</sup>	32.02 ± 2.85 <sup>ab</sup>	30.19 ± 2.77 <sup>b</sup>	<0.0001
(mg/dl)	0.81 ± 0.03 <sup>b</sup>	0.82 ± 0.03 <sup>b</sup>	$0.81 \pm 0.04^{b}$	$0.90 \pm 0.08^{b}$	0.84 ± 0.07 <sup>b</sup>	0.001
Systolic blood pressure (mmHg)	113.73 ± 3.11 <sup>b</sup>	119.17 ± 3.43ªb	109.13 ± 3.59 <sup>b</sup>	133.93 ± 2.75ª	135.00 ± 3.87ª	<0.0001
Diastolic blood pressure (mmHg)	71.79 ± 1.81 <sup>b</sup>	76.89 ± 1.58 <sup>ab</sup>	70.75 ± 3.25 <sup>b</sup>	85.19 ± 1.25ª	80.97 ± 2.23ª	<0.0001